

Collectin-11

Catalog No: #AP81825



Package Size: #AP81825-1 50ug #AP81825-2 100ug #AP81825-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collectin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Collectin-2,CL-2,CL2
Accession No.	Uniprot ID:Q2LK94
Uniprot	Q2LK94
GeneID	428649;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COLEC11

Note: This product is for in vitro research use only