

Collagen alpha-1(XXVII) chain

Catalog No: #AP81827



Package Size: #AP81827-1 50ug #AP81827-2 100ug #AP81827-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XXVII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KIAA1870
Accession No.	Uniprot ID:Q8IZC6
Uniprot	Q8IZC6
GeneID	85301;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: COL27A1

Note: This product is for in vitro research use only