

Claudin-1

Catalog No: #AP81828



Package Size: #AP81828-1 50ug #AP81828-2 100ug #AP81828-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Claudin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Senescence-associated epithelial membrane protein,CLD1,SEMP1,UNQ481/PRO944
Accession No.	Uniprot ID:O95832
Uniprot	O95832
GeneID	9076;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CLDN1

Note: This product is for in vitro research use only