

Cystatin-9

Catalog No: #AP81831

Package Size: #AP81831-1 50ug #AP81831-2 100ug #AP81831-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cystatin-like molecule,CLM,CTES7A
Accession No.	Uniprot ID:Q5W186
Uniprot	Q5W186
GeneID	128822;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CST9

Note: This product is for in vitro research use only