

C-C motif chemokine 20

Catalog No: #AP81834



Package Size: #AP81834-1 50ug #AP81834-2 100ug #AP81834-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta-chemokine exodus-1,CC chemokine LARC,CC chemokine ST38,Liver and activation-regulated chemokine,Macrophage inflammatory protein 3 alpha,MIP-3-alpha,Small-inducible cytokine A20,Larc,Scya20
Accession No.	Uniprot ID:O89093
Uniprot	O89093
GeneID	20297;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl20

Note: This product is for in vitro research use only