

Chitinase-3-like protein 1

Catalog No: #AP81839



Package Size: #AP81839-1 50ug #AP81839-2 100ug #AP81839-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chitinase-3-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	39 kDa synovial protein, Cartilage glycoprotein 39, CGP-39, YKL-40
Accession No.	Uniprot ID:P36222
Uniprot	P36222
GeneID	1116;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CHI3L1

Note: This product is for in vitro research use only