

Chromogranin-A

Catalog No: #AP81841



Package Size: #AP81841-1 50ug #AP81841-2 100ug #AP81841-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chromogranin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CgA,Pituitary secretory protein I,SP-I
Accession No.	Uniprot ID:P10645
Uniprot	P10645
GeneID	1113;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHGA

Note: This product is for in vitro research use only