

Chondroadherin

Catalog No: #AP81844



Package Size: #AP81844-1 50ug #AP81844-2 100ug #AP81844-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chondroadherin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cartilage leucine-rich protein
Accession No.	Uniprot ID:O55226
Uniprot	O55226
GeneID	12643;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Chad

Note: This product is for in vitro research use only