

JAK2 Mouse Monoclonal Antibody(3F10)

Catalog No: #44247

Package Size: #44247-1 50ul #44247-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | JAK2 Mouse Monoclonal Antibody(3F10) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 3F10 |
| Isotype | IgG1 |
| Applications | IHC |
| Species Reactivity | Hu Ms Rt |
| Specificity | The antibody detects endogenous JAK2 protein |
| Immunogen Type | Peptide |
| Immunogen Description | Synthetic Peptide of JAK2 |
| Target Name | JAK2 |
| Other Names | Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2) |
| Accession No. | Swiss-Prot:O60674Gene ID:3717 |
| Uniprot | O60674 |
| GeneID | 3717; |
| Calculated MW | 120kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only