

## Chromodomain-helicase-DNA-binding protein 8

Catalog No: #AP81851



Package Size: #AP81851-1 50ug #AP81851-2 100ug #AP81851-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chromodomain-helicase-DNA-binding protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CHD-8,ATP-dependent helicase CHD8,Axis duplication inhibitor,Duplin,Kiaa1564
Accession No.	Uniprot ID:Q09XV5
Uniprot	Q09XV5
GeneID	67772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Chd8

Note: This product is for in vitro research use only