

Collagen alpha-1(I) chain

Catalog No: #AP81854



Package Size: #AP81854-1 50ug #AP81854-2 100ug #AP81854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Collagen alpha-1(I) chain |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Alpha-1 type I collagen |
| Accession No. | Uniprot ID:Q9XSJ7 |
| Uniprot | Q9XSJ7 |
| GeneID | 403651; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: COL1A1

Note: This product is for in vitro research use only