

Adapter molecule crk

Catalog No: #AP81855



Package Size: #AP81855-1 50ug #AP81855-2 100ug #AP81855-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adapter molecule crk
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Proto-oncogene c-Crk,p38
Accession No.	Uniprot ID:Q04929
Uniprot	Q04929
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CRK

Note: This product is for in vitro research use only