

## Versican core protein

Catalog No: #AP81857



Package Size: #AP81857-1 50ug #AP81857-2 100ug #AP81857-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Versican core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chondroitin sulfate proteoglycan core protein 2,Chondroitin sulfate proteoglycan 2,Glial hyaluronate-binding protein,GHAP,Large fibroblast proteoglycan,PG-M,Cspg2
Accession No.	Uniprot ID:Q9ERB4
Uniprot	Q9ERB4
GeneID	114122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Vcan

Note: This product is for in vitro research use only