

## JAK2 Mouse Monoclonal Antibody(9B4)

Catalog No: #44248

Package Size: #44248-1 50ul #44248-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	JAK2 Mouse Monoclonal Antibody(9B4)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9B4
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous JAK2 protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of JAK2
Target Name	JAK2
Other Names	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Accession No.	Swiss-Prot:O60674 Gene ID:3717
Uniprot	O60674
GeneID	3717;
Calculated MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK2 Mouse mAb diluted at 1:200

---

Note: This product is for in vitro research use only