

H(+)/Cl(-) exchange transporter 5

Catalog No: #AP81861

Package Size: #AP81861-1 50ug #AP81861-2 100ug #AP81861-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	H(+)/Cl(-) exchange transporter 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Chloride channel protein 5, CIC-5, Chloride transporter CIC-5
Accession No.	Uniprot ID:Q9GKE7
Uniprot	Q9GKE7
GeneID	397263;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CLCN5

Note: This product is for in vitro research use only