

T-cell surface glycoprotein CD3 epsilon chain

Catalog No: #AP81863



Package Size: #AP81863-1 50ug #AP81863-2 100ug #AP81863-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD3 epsilon chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CD3e
Accession No.	Uniprot ID:P27597
Uniprot	P27597
GeneID	442981;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD3E

Note: This product is for in vitro research use only