

T-cell surface glycoprotein CD8 alpha chain

Catalog No: #AP81864



Package Size: #AP81864-1 50ug #AP81864-2 100ug #AP81864-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD8 alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD8 antigen 32 kDa chain, OX-8 membrane antigen, CD8a
Accession No.	Uniprot ID:P07725
Uniprot	P07725
GeneID	24930;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd8a

Note: This product is for in vitro research use only