

C-X-C motif chemokine 2

Catalog No: #AP81869



Package Size: #AP81869-1 50ug #AP81869-2 100ug #AP81869-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokine-induced neutrophil chemoattractant 3,CINC-3,Macrophage inflammatory protein 2,MIP2,Cinc3,Mip-2,Mip2,Scyb2
Accession No.	Uniprot ID:P30348
Uniprot	P30348
GeneID	114105;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcl2

Note: This product is for in vitro research use only