

C-X-C motif chemokine 9

Catalog No: #AP81870



Package Size: #AP81870-1 50ug #AP81870-2 100ug #AP81870-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gamma-interferon-induced monokine, Monokine induced by interferon-gamma, HuMIG, Small-inducible cytokine B9, CMK, MIG, SCYB9
Accession No.	Uniprot ID:Q07325
Uniprot	Q07325
GeneID	4283;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CXCL9

Note: This product is for in vitro research use only