

Cytochrome P450 2E1

Catalog No: #AP81878



Package Size: #AP81878-1 50ug #AP81878-2 100ug #AP81878-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2E1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	4-nitrophenol 2-hydroxylase,CYPIIE1,CYP2E
Accession No.	Uniprot ID:O18963
Uniprot	O18963
GeneID	282213;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP2E1

Note: This product is for in vitro research use only