

CMRF35-like molecule 7

Catalog No: #AP81883



Package Size: #AP81883-1 50ug #AP81883-2 100ug #AP81883-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CMRF35-like molecule 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CLM-7,CD300 antigen-like family member B,Immune receptor expressed on myeloid cells 3,IREM-3,Leukocyte mono-Ig-like receptor 5,CD300b,Irem3,Lmir5
Accession No.	Uniprot ID:Q3U497
Uniprot	Q3U497
GeneID	217304;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd300lb

Note: This product is for in vitro research use only