

Cytochrome P450 2D4

Catalog No: #AP81887



Package Size: #AP81887-1 50ug #AP81887-2 100ug #AP81887-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2D4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CYP11D4,Cytochrome P450-CMF3,Cytochrome P450-DB4,Debrisoquine 4-hydroxylase,Cyp2d-4
Accession No.	Uniprot ID:P13108
Uniprot	P13108
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp2d4

Note: This product is for in vitro research use only