

C-type lectin domain family 4 member E

Catalog No: #AP81891



Package Size: #AP81891-1 50ug #AP81891-2 100ug #AP81891-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-type lectin domain family 4 member E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-type lectin superfamily member 9,Macrophage-inducible C-type lectin,Clecsf9,Mincle
Accession No.	Uniprot ID:Q67EQ1
Uniprot	Q67EQ1
GeneID	450223;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Clec4e

Note: This product is for in vitro research use only