

Complement component C7

Catalog No: #AP81892

Package Size: #AP81892-1 50ug #AP81892-2 100ug #AP81892-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Complement component C7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q29RQ1
Uniprot	Q29RQ1
GeneID	507339;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C7

Note: This product is for in vitro research use only