

Collagen alpha-1(XVII) chain

Catalog No: #AP81894



Package Size: #AP81894-1 50ug #AP81894-2 100ug #AP81894-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XVII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	180 kDa bullous pemphigoid antigen 2,Bullous pemphigoid antigen 2,BP180,BPAG2
Accession No.	Uniprot ID:Q90584
Uniprot	Q90584
GeneID	396503;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL17A1

Note: This product is for in vitro research use only