

H(+)/Cl(-) exchange transporter 5

Catalog No: #AP81896

Package Size: #AP81896-1 50ug #AP81896-2 100ug #AP81896-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	H(+)/Cl(-) exchange transporter 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Chloride channel protein 5,ClC-5,Chloride transporter ClC-5,Clc5
Accession No.	Uniprot ID:Q9WVD4
Uniprot	Q9WVD4
GeneID	12728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Clcn5

Note: This product is for in vitro research use only