

Calponin-1

Catalog No: #AP81897



Package Size: #AP81897-1 50ug #AP81897-2 100ug #AP81897-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calponin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Basic calponin, Calponin H1, smooth muscle
Accession No.	Uniprot ID:Q2HJ38
Uniprot	Q2HJ38
GeneID	534583;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CNN1

Note: This product is for in vitro research use only