

Complement factor D

Catalog No: #AP81902



Package Size: #AP81902-1 50ug #AP81902-2 100ug #AP81902-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	28 kDa adipocyte protein,Adipsin,C3 convertase activator,Properdin factor D,Adn,Df
Accession No.	Uniprot ID:P03953
Uniprot	P03953
GeneID	11537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cfd

Note: This product is for in vitro research use only