

Granulocyte colony-stimulating factor receptor

Catalog No: #AP81908



Package Size: #AP81908-1 50ug #AP81908-2 100ug #AP81908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Granulocyte colony-stimulating factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	G-CSF receptor,CD114,Csfgr
Accession No.	Uniprot ID:P40223
Uniprot	P40223
GeneID	12986;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Csf3r

Note: This product is for in vitro research use only