

Complement C4-A

Catalog No: #AP81910



Package Size: #AP81910-1 50ug #AP81910-2 100ug #AP81910-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement C4-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acidic complement C4, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2, CO4, CPAMD2
Accession No.	Uniprot ID:P0C0L4
Uniprot	P0C0L4
GeneID	720;721;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C4A

Note: This product is for in vitro research use only