

Cytochrome P450 2C9

Catalog No: #AP81913



Package Size: #AP81913-1 50ug #AP81913-2 100ug #AP81913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Cytochrome P450 2C9 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | (R)-limonene 6-monooxygenase,(S)-limonene 6-monooxygenase,(S)-limonene 7-monooxygenase,CYP11C9,Cholesterol 25-hydroxylase,Cytochrome P-450MP,Cytochrome P450 MP-4,Cytochrome P450 MP-8,Cytochrome P450 PB-1,S-mephenytoin 4-hydroxylase,CYP2C10 |
| Accession No. | Uniprot ID:P11712 |
| Uniprot | P11712 |
| GeneID | 1559; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP2C9

Note: This product is for in vitro research use only