

Metal transporter CNNM2

Catalog No: #AP81914



Package Size: #AP81914-1 50ug #AP81914-2 100ug #AP81914-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metal transporter CNNM2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ancient conserved domain-containing protein 2,Cyclin-M2,Acdp2
Accession No.	Uniprot ID:Q5U2P1
Uniprot	Q5U2P1
GeneID	294014;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cnnm2

Note: This product is for in vitro research use only