

C-C chemokine receptor type 7

Catalog No: #AP81919



Package Size: #AP81919-1 50ug #AP81919-2 100ug #AP81919-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-C CKR-7,BLR2,CDw197,Epstein-Barr virus-induced G-protein coupled receptor 1,EBI1,MIP-3 beta receptor,CD197,CMKBR7,EBI1,EVI1
Accession No.	Uniprot ID:P32248
Uniprot	P32248
GeneID	1236;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCR7

Note: This product is for in vitro research use only