

Catenin delta-1

Catalog No: #AP81920

Package Size: #AP81920-1 50ug #AP81920-2 100ug #AP81920-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Catenin delta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cadherin-associated Src substrate,CAS,p120 catenin,p120(ctn),p120(cas),Catns,Kiaa0384
Accession No.	Uniprot ID:P30999
Uniprot	P30999
GeneID	12388;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ctnnd1

Note: This product is for in vitro research use only