

Cullin-5

Catalog No: #AP81921



Package Size: #AP81921-1 50ug #AP81921-2 100ug #AP81921-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cullin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CUL-5,Vasopressin-activated calcium-mobilizing receptor 1,VACM-1,Vacm1
Accession No.	Uniprot ID:Q9JJ31
Uniprot	Q9JJ31
GeneID	64624;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cul5

Note: This product is for in vitro research use only