

## C-X-C motif chemokine 11

Catalog No: #AP81922



Package Size: #AP81922-1 50ug #AP81922-2 100ug #AP81922-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-X-C motif chemokine 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interferon-inducible T-cell alpha chemoattractant,I-TAC,Small-inducible cytokine B11,Scyb11
Accession No.	Uniprot ID:Q9JHH5
Uniprot	Q9JHH5
GeneID	56066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cxcl11

Note: This product is for in vitro research use only