

2',3'-cyclic-nucleotide 3'-phosphodiesterase

Catalog No: #AP81924



Package Size: #AP81924-1 50ug #AP81924-2 100ug #AP81924-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | 2',3'-cyclic-nucleotide 3'-phosphodiesterase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | CNP,Cnp1 |
| Accession No. | Uniprot ID:P16330 |
| Uniprot | P16330 |
| GeneID | 12799; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cnp

Note: This product is for in vitro research use only