

Cannabinoid receptor 2

Catalog No: #AP81926



Package Size: #AP81926-1 50ug #AP81926-2 100ug #AP81926-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cannabinoid receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CB-2
Accession No.	Uniprot ID:P47936
Uniprot	P47936
GeneID	12802;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cnr2

Note: This product is for in vitro research use only