

Ciliary neurotrophic factor

Catalog No: #AP81927



Package Size: #AP81927-1 50ug #AP81927-2 100ug #AP81927-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ciliary neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CNTF
Accession No.	Uniprot ID:P51642
Uniprot	P51642
GeneID	12803;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cntf

Note: This product is for in vitro research use only