

Cytochrome P450 19A1

Catalog No: #AP81935



Package Size: #AP81935-1 50ug #AP81935-2 100ug #AP81935-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 19A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Aromatase,CYPXIX,Cytochrome P-450AROM,Estrogen synthase,CYP19
Accession No.	Uniprot ID:Q5QQX7
Uniprot	Q5QQX7
GeneID	494003;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP19A1

Note: This product is for in vitro research use only