

Complement receptor type 2

Catalog No: #AP81936



Package Size: #AP81936-1 50ug #AP81936-2 100ug #AP81936-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cr2, Complement C3d receptor, Epstein-Barr virus receptor, EBV receptor, CD21, C3DR
Accession No.	Uniprot ID:P20023
Uniprot	P20023
GeneID	1380;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CR2

Note: This product is for in vitro research use only