

## Calcium release-activated calcium channel protein 1

Catalog No: #AP81937



Package Size: #AP81937-1 50ug #AP81937-2 100ug #AP81937-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Calcium release-activated calcium channel protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Protein orai-1, Transmembrane protein 142A, CRACM1, TMEM142A, RCJMB04_8f15
Accession No.	Uniprot ID:Q5ZL05
Uniprot	Q5ZL05
GeneID	416845;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ORAI1

Note: This product is for in vitro research use only