

Transcriptional repressor CTCFL

Catalog No: #AP81941



Package Size: #AP81941-1 50ug #AP81941-2 100ug #AP81941-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcriptional repressor CTCFL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brother of the regulator of imprinted sites,CCCTC-binding factor,CTCF paralog,CTCF-like protein,Boris
Accession No.	Uniprot ID:A2APF3
Uniprot	A2APF3
GeneID	664799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ctcf

Note: This product is for in vitro research use only