

Choline transporter-like protein 4

Catalog No: #AP81942



Package Size: #AP81942-1 50ug #AP81942-2 100ug #AP81942-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Choline transporter-like protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Solute carrier family 44 member 4,Ctl4,Ng22
Accession No.	Uniprot ID:Q6MG71
Uniprot	Q6MG71
GeneID	294255;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc44a4

Note: This product is for in vitro research use only