

Cytosolic non-specific dipeptidase

Catalog No: #AP81943



Package Size: #AP81943-1 50ug #AP81943-2 100ug #AP81943-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic non-specific dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CNDP dipeptidase 2
Accession No.	Uniprot ID:Q6Q0N1
Uniprot	Q6Q0N1
GeneID	291394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cndp2

Note: This product is for in vitro research use only