

## Metal transporter CNNM2

Catalog No: #AP81945



Package Size: #AP81945-1 50ug #AP81945-2 100ug #AP81945-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metal transporter CNNM2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ancient conserved domain-containing protein 2,mACDP2,Cyclin-M2,Acdp2
Accession No.	Uniprot ID:Q3TWN3
Uniprot	Q3TWN3
GeneID	94219;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cnnm2

Note: This product is for in vitro research use only