

Carcinoembryonic antigen-related cell adhesion molecule 1

Catalog No: #AP81950

Package Size: #AP81950-1 50ug #AP81950-2 100ug #AP81950-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carcinoembryonic antigen-related cell adhesion molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Biliary glycoprotein 1,BGP-1,Biliary glycoprotein D,MHVR1,Murine hepatitis virus receptor,MHV-R,CD66a,Bgp,Bgp1
Accession No.	Uniprot ID:P31809
Uniprot	P31809
GeneID	26365;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ceacam1

Note: This product is for in vitro research use only