

# 10 kDa heat shock protein, mitochondrial

Catalog No: #AP81953



Package Size: #AP81953-1 50ug #AP81953-2 100ug #AP81953-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	10 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hsp10,10 kDa chaperonin,Chaperonin 10,CPN10,Early-pregnancy factor,EPF
Accession No.	Uniprot ID:P61604
Uniprot	P61604
GeneID	3336;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HSPE1

Note: This product is for in vitro research use only