

Cytochrome P450 2B6

Catalog No: #AP81957



Package Size: #AP81957-1 50ug #AP81957-2 100ug #AP81957-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2B6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	1,4-cineole 2-exo-monooxygenase,CYP11B6,Cytochrome P450 I1B1
Accession No.	Uniprot ID:P20813
Uniprot	P20813
GeneID	1555;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP2B6

Note: This product is for in vitro research use only