

Metal transporter CNNM2

Catalog No: #AP81963



Package Size: #AP81963-1 50ug #AP81963-2 100ug #AP81963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metal transporter CNNM2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ancient conserved domain-containing protein 2,Cyclin-M2,ACDP2
Accession No.	Uniprot ID:Q9H8M5
Uniprot	Q9H8M5
GeneID	54805;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CNNM2

Note: This product is for in vitro research use only